

MECHANICAL ULTRASONIC DESKTOP CLEANER ULC14-32



MECHANICAL ULTRASONIC DESKTOP CLEANER

ULC14-32

Mechanical Ultrasonic desktop cleaner has high capacity for internal operations, high cleaning accuracy Effectively cleans small dirt particles, even deep holes, small seams and concealment of workpieces.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC14-32 MECHANICAL ULTRASONIC DESKTOP CLEANER

Mechanical temperature and time setting

Industrial high Q transducers, cleaning conversion efficiency is higher

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control



SPECIFICATIONS

Model	ULC14-32
Capacity	3.2 L
Working frequency	40 KHz
Ultrasonic Power	120 W
Heating power	100 W
Timing time	1- 60 mins
Temperature range	0- 80°C
Cover	Yes
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	240x135x100 mm
Packing size	342x240x310 mm
Gross weight	4.77 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com