

# HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER ULC15-11



# HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER ULC15-11

High frequency ultrasonic cleaner is mainly used in medical industry, computers, microelectronics computers. It uses high frequency to clean the sub micron dirt without damaging the equipment

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC15-11 HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER

Digital display and control ultrasonic power adjustable range from 0 to 100%

Special high frequency ultrasonic transducer (53/68/80/100/132 KHz)

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control



### SPECIFICATIONS

Model	ULC15-11
Capacity	11 L
Working frequency	53/68/80/100/132 KHz
Ultrasonic Power	240 W
Heating power	300 W
Timing time	1- 99 mins
Temperature range	0- 80°C
Power regulation	0-100 %
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	300x240x150 mm
Packing size	300x240x150 mm
Gross weight	8.5 kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)