

HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER ULC15-22



HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER ULC15-22

High frequency ultrasonic cleaner is mainly used in medical industry, computers, microelectronics computers. It uses high frequency to clean the sub micron dirt without damaging the equipment

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC15-22 HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER

Digital display and control ultrasonic power adjustable range from 0 to 100%

Special high frequency ultrasonic transducer (53/68/80/100/132 Khz)

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control



SPECIFICATIONS

Model	ULC15-22
Capacity	22 L
Working frequency	53/68/80/100/132 Khz
Ultrasonic Power	480 W
Heating power	600 W/220V 500 W/110V
Timing time	1- 99 mins
Temperature range	0- 80°C
Power regulation	0-100 %
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	500x300x150 mm
Packing size	500x300x150 mm
Gross weight	14.51 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com