

# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER ULC17-10A



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

## ULC17-10A

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC17-10A SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz



### SPECIFICATIONS

Model	ULC17-10A
Capacity	10 L
Ultrasonic frequency	80 KHz
Time range/ Heating temp	1-480 mins / RT-80°C
Ultrasonic power	300 W
Heating power	400 W
Tank size	300x240x150 mm
Packing Size	446x346x463 mm
Gross Weight	12.1 Kg

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)