

# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER ULC17-15



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

## ULC17-15

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC17-15 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz



## SPECIFICATIONS

|                          |                      |
|--------------------------|----------------------|
| Model                    | ULC17-15             |
| Capacity                 | 15 L                 |
| Ultrasonic frequency     | 40 Khz               |
| Time range/ Heating temp | 1-30 mins / 20- 80°C |
| Ultrasonic power         | 360 W                |
| Heating power            | 600 W                |
| Tank size                | 330x300x150 mm       |
| Packing Size             | 460x430x390 mm       |
| Gross Weight             | 10 Kg                |

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)