

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER ULC17-15



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

ULC17-15

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC17-15 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz



SPECIFICATIONS

Model	ULC17-15
Capacity	15 L
Ultrasonic frequency	40 Khz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	360 W
Heating power	600 W
Tank size	330x300x150 mm
Packing Size	460x430x390 mm
Gross Weight	10 Kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com