

# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER ULC17-27F



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

## ULC17-27F

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC17-27F SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz



### SPECIFICATIONS

Model	ULC17-27F
Capacity	27 L
Ultrasonic frequency	80 KHz
Time range/ Heating temp	1-480 mins /0- 80°C
Ultrasonic power	700 W
Heating power	800 W
Tank size	500x300x180 mm
Packing Size	1100x690x820 mm
Gross Weight	101 Kg

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)