

# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER ULC17-30



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

## ULC17-30

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC17-30 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz



## SPECIFICATIONS

Model	ULC17-30
Capacity	30 L
Ultrasonic frequency	40 Khz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	600 W
Heating power	1000 W
Tank size	500x300x200 mm
Packing Size	630x430x430 mm
Gross Weight	15 Kg

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)