

DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER ULC18-6A



DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER ULC18-6A

Double adjustable frequencies ultrasonic cleaner uses frequency of 45 KHz /80 KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC18-6A DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

Over- Voltage and over current protection

Standard basket and lid

Low level water protection

Fully microprocessor controller and programmable

Time and temperature real time display

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz



SPECIFICATIONS

| | |
|---------------------------------|----------------------|
| Model | ULC18-6A |
| Capacity | 6 L |
| Adjustable ultrasonic frequency | 45 KHz /80 KHz |
| Time range/ Heating temp | 1-480 mins / RT-80°C |
| Ultrasonic power | 200 W |
| Heating power | 400 W |
| Tank size | 300x150x150 mm |
| Packing Size | 445x250x430 mm |
| Gross Weight | 9.5 Kg |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com