

# ULTRASONIC CLEANER ULC19-10C



# ULTRASONIC CLEANER ULC19-10C

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC19-10C ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



## SPECIFICATIONS

Model	ULC19-10C
Capacity	10 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	240 W
Heating power	500 W
Tank size	300x240x150 mm
External Size	300x270x310 mm

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)