

LABSTAC

ULTRASONIC CLEANER ULC19-15A



labstac.com

ULTRASONIC CLEANER ULC19-15A

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC19-15A ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



SPECIFICATIONS

Model	ULC19-15A
Capacity	15 L
Ultrasonic frequency	Single frequency, 40 Khz
Time range/ Heating temp	1-60 min/RT-80°C
Ultrasonic power	360 W
Heating power	500 W
Tank size	330x300x150 mm
External Size	360x330x310 mm
Packing Size	460x400x400 mm
Gross weight	10 Kg
Power Supply	AC 220V, 50Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com