

# ULTRASONIC CLEANER ULC19-15B



# ULTRASONIC CLEANER ULC19-15B

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC19-15B ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



## SPECIFICATIONS

Model	ULC19-15B
Capacity	15 L
Ultrasonic frequency	Double frequency 28/40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	360 W
Heating power	500 W
Tank size	330x300x150 mm
External Size	360x330x310 mm

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)