

# ULTRASONIC CLEANER ULC19-225A



# ULTRASONIC CLEANER ULC19-225A

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC19-225A ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



## SPECIFICATIONS

Model	ULC19-225A
Capacity	22.5 L
Ultrasonic frequency	Single frequency, 40 Khz
Time range/ Heating temp	1-60 min/RT-80°C
Ultrasonic power	480 W
Heating power	500 W
Tank size	500x300x150 mm
External Size	530x330x310 mm
Packing Size	640x400x400 mm
Gross weight	13 Kg
Power Supply	AC 220V, 50Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)