

LABSTAC

ULTRASONIC CLEANER ULC19-225C



labstac.com

ULTRASONIC CLEANER ULC19-225C

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC19-225C ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



SPECIFICATIONS

Model	ULC19-225C
Capacity	22.5 L
Ultrasonic frequency	Single frequency, 40 Khz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	480 W
Heating power	500 W
Tank size	500x300x150 mm
External Size	530x330x310 mm

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com