

# ULTRASONIC CLEANER ULC19-225C



# ULTRASONIC CLEANER ULC19-225C

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC19-225C ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



## SPECIFICATIONS

|                          |                          |
|--------------------------|--------------------------|
| Model                    | ULC19-225C               |
| Capacity                 | 22.5 L                   |
| Ultrasonic frequency     | Single frequency, 40 Khz |
| Time range/ Heating temp | 1-99 min/RT-80°C         |
| Ultrasonic power         | 480 W                    |
| Heating power            | 500 W                    |
| Tank size                | 500x300x150 mm           |
| External Size            | 530x330x310 mm           |

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)