

# ULTRASONIC CLEANER ULC19-30C



# ULTRASONIC CLEANER ULC19-30C

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC19-30C ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



## SPECIFICATIONS

Model	ULC19-30C
Capacity	30 L
Ultrasonic frequency	Single frequency, 40 Khz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	600 W
Heating power	500 W
Tank size	500x300x200 mm
External Size	530x330x350 mm

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)