

ULTRASONIC CLEANER ULC19-30C



ULTRASONIC CLEANER ULC19-30C

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC19-30C ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



SPECIFICATIONS

| | |
|--------------------------|--------------------------|
| Model | ULC19-30C |
| Capacity | 30 L |
| Ultrasonic frequency | Single frequency, 40 KHz |
| Time range/ Heating temp | 1-99 min/RT-80°C |
| Ultrasonic power | 600 W |
| Heating power | 500 W |
| Tank size | 500x300x200 mm |
| External Size | 530x330x350 mm |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com