

ULTRASONIC CLEANER ULC19-65



ULTRASONIC CLEANER ULC19-65

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC19-65 ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



SPECIFICATIONS

Model	ULC19-65
Capacity	6.5 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	180 W
Heating power	500 W
Tank size	300x150x150 mm
External Size	330x180x310 mm
Packing Size	430x250x350 mm
Gross weight	6 Kg
Power Supply	AC 220V, 50Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com