

LABSTAC

ULTRASONIC CLEANER ULC19-65C



labstac.com

ULTRASONIC CLEANER ULC19-65C

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC19-65C ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



SPECIFICATIONS

Model	ULC19-65C
Capacity	6.5 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	180 W
Heating power	500 W
Tank size	300x150x150 mm
External Size	300x180x310 mm

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com