

**-10°C TO -25°C ULTRA
LOW TEMPERATURE
UPRIGHT FREEZER
ULT13-100**



-10°C TO -25°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER ULT13-100

Engineered for storage of number of samples like test samples and plasma to biologicals, enzymes for genetic research, sensitive materials and temperature. Features like low energy consumption, low noise level, low heat dissipation and safety towards audible alarm for high/ low temperature makes this product perfect.

Used in Laboratories, Routine Cold Storage, Biologicals, Pharmaceuticals, Clinical, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

ULT13-100 -10°C TO -25°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER

Undercounter Mini type Drawer x 4



SPECIFICATIONS

Model	ULT13-100
Capacity	100 L
Temperature Range	-10°C to -25°C
Alarm	High&low temp, power failure, control error, door open
External Dimension	540X540X880 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com