

# -10°C TO -25°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER ULT13-358



# -10°C TO -25°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER ULT13-358

Engineered for storage of number of samples like test samples and plasma to biologicals, enzymes for genetic research, sensitive materials and temperature. Features like low energy consumption, low noise level, low heat dissipation and safety towards audible alarm for high/ low temperature makes this product perfect.

Used in Laboratories, Routine Cold Storage, Biologicals, Pharmaceuticals, Clinical, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

## ULT13-358 -10°C TO -25°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER

Double Door Drawer x 12 USB port



### SPECIFICATIONS

Model	ULT13-358
Capacity	358 L
Temperature Range	-10°C to -25°C
Alarm	High&low temp, power failure, control error, door open, high ambient temperture
External Dimension	720×830×1870 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)