

BIOSAFETY TRANSPORT BOX ULT92-30




BIOSAFETY TRANSPORT BOX ULT92-30

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc. Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

ULT92-30 BIOSAFETY TRANSPORT BOX

AC & DC available Able to used in cars



SPECIFICATIONS

Model	ULT92-30
Capacity	30 L
Temperature Range	-25 to 30°C
Alarm	High&low temp, control error, power failure
External Dimension	610x360x460 mm
Power	DC12V; AC100~220V, 50/60Hz
Optional Parts	Basket; Datalogger

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com