

# BIOSAFETY TRANSPORT BOX ULT92-80



# BIOSAFETY TRANSPORT BOX ULT92-80

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc. Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

## ULT92-80 BIOSAFETY TRANSPORT BOX



### SPECIFICATIONS

|                    |   |
|--------------------|---|
| Model              | ULT92-80                                    |
| Capacity           | 80 L  |
| Temperature Range  | -25 to 30°C                                 |
| Alarm              | High&low temp, control error, power failure |
| External Dimension | 930x580x570 mm                              |
| Power              | DC12V; AC100~220V, 50/60Hz                  |
| Optional Parts     | Basket; Datalogger                          |

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)