

# PORTABLE AUTOCLAVE



# PORTABLE AUTOCLAVE

Portable Autoclave is used for general purpose sterilization. Saturated steam is used to sterilize articles rapidly and efficiently. It is also used as a cooker in plateau areas and also for making high quality drinking water in Drinking Water Plant.

Used in Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research..

Also known as Portable Sterilizer, Steam Autoclave Sterilizer, Portable Autoclave, Pressure Steam Sterilizer, Portable Pressure Steam Sterilizer, Dental Autoclave.

## AUT41 PORTABLE AUTOCLAVE



### SPECIFICATIONS

| Model                     | AUT41-18                | AUT41-24       |
|---------------------------|-------------------------|----------------|
| Capacity                  | 18 L                    | 24 L           |
| Function                  | Electric or LPG heated. |                |
| Sterilization Temperature | 126°C                   |                |
| Working Pressure          | 0.14-0.16 Mpa           |                |
| Maximum Working Pressure  | 0.165 Mpa               |                |
| Chamber Dimension         | 280x260 mm              | 280x390mm      |
| Overall Dimension (LxWxH) | 410x410x430 mm          | 410x410x550 mm |
| Chamber Material          | Stainless Steel         |                |
| Weight                    | 14/16 kg                | 15/17 kg       |
| Power                     | 2 KW                    |                |
| Power Supply              | 220 V, 50 Hz            |                |

For Models AUT41-18

For Models AUT41-24



AUT41-18



AUT41-24

## AUT42 PORTABLE AUTOCLAVE



### SPECIFICATIONS

| Model                        | AUT42-18        | AUT42-24       |
|------------------------------|-----------------|----------------|
| Capacity                     | 18 L            | 24 L           |
| Function                     | Electric heated |                |
| Sterilization Temperature    | 126°C           |                |
| Steam Exhausting Temperature | ≤1°C            |                |
| Working Pressure             | 0.14-0.16Mpa    |                |
| Maximum Working Pressure     | 0.165 Mpa       |                |
| Chamber Dimension            | 280x260 mm      | 280x390mm      |
| Overall Dimension (LxWxH)    | 410x410x520 mm  | 410x410x630 mm |
| Chamber Material             | Stainless Steel |                |
| Weight                       | 14/16 kg        | 15/17 kg       |
| Power                        | 2 KW            |                |
| Power Supply                 | 220 V, 50 Hz    |                |

For Models AUT42-18

For Models AUT42-24



AUT42-18



AUT42-24

## AUT43 PORTABLE AUTOCLAVE

### SPECIFICATIONS

| Model                        | AUT43-18        | AUT43-24       |
|------------------------------|-----------------|----------------|
| Capacity                     | 18 L            | 24 L           |
| Function                     | Electric heated |                |
| Sterilization Temperature    | 126°C           |                |
| Steam Exhausting Temperature | ≤1°C            |                |
| Working Pressure             | 0.14-0.16Mpa    |                |
| Maximum Working Pressure     | 0.165 Mpa       |                |
| Timer Range                  | 0-99min         |                |
| Chamber Dimension            | 280x260mm       | 270x500mm      |
| Overall Dimension (LxWxH)    | 440x420x630 mm  | 440x420x750 mm |
| Chamber Material             | Stainless Steel |                |
| Weight                       | 12/14 kg        | 14/15 kg       |
| Power                        | 2 KW            |                |
| Power Supply                 | 220 V, 50 Hz    |                |

For Models AUT43-18

For Models AUT43-24

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)