

HIGH TEMPERATURE CIRCULATOR



HIGH TEMPERATURE CIRCULATOR

BAC037 HIGH TEMPERATURE CIRCULATOR

Microprocessor control with timing function
Digital display
Audible and visible alarm for over temperature



SPECIFICATIONS

Model	BAC037-22	BAC037-30
Capacity	22 L	30 L
Electrical Requirement	220 V 50 HZ	
Power Consumption	1000 W	
Temperature Range	RT+5-70°C	
Temperature Stability	±0.3°C	
Display Resolution	0.1 °C	
Interior dimension	500Wx300Dx150H mm	600Wx300Dx190H mm
Exterior dimension	660Wx340Dx200H mm	770Wx370Dx280H mm
Timing Range	1~999 min	



BAC037-22



BAC037-30

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com