

AIR SHOWER PASS BOX



AIR SHOWER PASS BOX

BSC51 AIR SHOWER BOX

Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness



SPECIFICATIONS

Model	BSC51-500	BSC51-600	BSC51-700
Clean Air Speed	0.36~0.54m/s		
UV Lamp	15Wx1		
External Size (WxDxH)	660x570x1000 mm	760x670x1100 mm	860x770x1200 mm
Internal Size (WxDxH)	500x500x500 mm	600x600x600 mm	700x700x700 mm
Material	304 stainless steel		
Interlock	Electronical		
HEPA Filter	Efficiency $\geq 99.995\%$ at $0.3\mu\text{m}$		
Gross Weight	80 kg	100 kg	120 kg
Packing Size(WxDxH)	780x700x1250 mm	880x800x1350mm	980x900x1450mm
Power Supply	220V, 50Hz; 110V, 50/60Hz		

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com