

PASS BOX



PASS BOX

BSC61 PASS BOX

Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness



SPECIFICATIONS

Model	BSC61-500	BSC61-600	BSC61-700
UV Lamp	8Wx1		
External Size (WxDxH)	660x570x580 mm	760x670x680 mm	860x770x780 mm
Internal Size (WxDxH)	500x500x500 mm	600x600x600 mm	700x700x700 mm
Material	304 stainless steel		
Interlock	Electronical		
Gross Weight	55 kg	70kg	90kg
Packing Size(WxDxH)	780x700x770 mm	880x800x870 mm	980x900x970 mm
Power Supply	220V, 50Hz; 110V, 50/60Hz		

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com