

BIOCHEMICAL INCUBATOR



BIOCHEMICAL INCUBATOR

Glass door for the easy visualization and compact with exclusive incubating performance in valuable tests. Made up of stainless steel which makes it durable. Biochemical Incubator is applicable in research institutes, university laboratories, environmental protection, forestry and animal husbandry industries for culture preservation, bacteria, and microbes. Used in Environmental protection, Forestry, Research institutes, University laboratories, Animal husbandry industries, Culture reservation, bacteria, and microbes.

Also known as Laboratory Biochemical Incubator.

INC41 BIOCHEMICAL INCUBATOR

ALLCOLD Refrigeration Technology: Auto defrost, multiple security system, long time running, environmental protection, high efficiency and energy saving, Auto-defrost function, Imported DuPont SUVA R134 a environmental refrigerant.

ALLFLOW Perfect Air Current Cycling: Perfect forced convection, maximum number of working room, minimum temperature recovery time after the opening, world famous axial fan, perfect air current design.

ALLSENS Programmable PID Control: Adaptive PID controller precisely controls the temperature and humidity, prevent temperature soaring, keep working room temperature stable and uniformity.

Excellent Imported temperature and humidity Sensor.

Efficient isolation Design



SPECIFICATIONS

| Model | INC41-100 | INC41-150 | INC41-250 | INC41-400 |
|----------------------------|------------------------------------|-------------------|--------------------|--------------------|
| Capacity | 100 L | 150 L | 250 L | 400 L |
| Temperature Range | 0-65°C | | | |
| Ambient Temperature | 10~30°C | | | |
| Temperature Fluctuation | ±0.5°C | | | |
| Temperature Uniformity | ±1 | | | ±1.5 |
| Temperature Resolution | 0.1°C | | | |
| Ambient Humidity | <70% | | | |
| Interior | SUS304 Stainless Steel Materials | - | - | - |
| Timer Range | 0-99 hrs, 0-9999 min, can be chose | | | |
| Convection Mode | Forced Convection | | | |
| Internal Dimension | 490Wx390Dx610H mm | 510Wx390Dx760H mm | 510Wx450Dx1100H mm | 600Wx640Dx1050H mm |
| Exterior Dimension (WxDxH) | 630x680x1250 mm | 650x680x1410 mm | 650x740x1730 mm | 745x930x1700 mm |
| Package Size | 1450x730x830 mm | 1605x760x830 mm | 1935x760x890 mm | 1907x815x1090 mm |
| Shelves/Trays | 3 | | 4 | |
| Weight | 95/120 kg | 107/133 kg | 135/162.5 kg | 158/158 kg |
| Power | 850 W | 1000 W | 1100 W | |
| Power Supply | Single phase AC220 V/50 Hz | | | |



INC41-100



INC41-150



INC41-250



INC41-400

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com