

# STRECHER



# STRECHER

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## HBD8A-PT PATIENT STRETCHER



### SPECIFICATIONS

Model	HBD8A-PT
Dimension	1860x590x700 mm

## HBD8C-00 HOSPITAL TRANSFER STRETCHER TROLLEY

It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;

Aluminium-Alloy protective rails, can be folded and hidden under the surface on no use;

Central-lock castors, diameter150mm. Matched with direction castor, help the operator controlling the direction easily;

Rotate the crand and handle on bedboard, patient on the mattress can be moved parallelly to another bed or flat table



### SPECIFICATIONS

Model	HBD8C-00
Height adjustment	920-610 mm
Net Weight	89 kg

## HBD8F-MS MANUAL TRANSFER STRETCHER

Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easy control the stretcher direction

Standard accessories: Mattress, I.V.Pole



### SPECIFICATIONS

Model	HBD8F-MS
Back-rest	0 – 85°
Height range (bed board to ground)	862 – 566 mm
Standard accessory	Mattress, IV pole
Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	2013x700 mm
Net Weight	82 kg

## HBD8G-MS1 MANUAL TRANSFER STRETCHER

Single crank for height adjustment Turn by clockwise, stretcher will be up; Turn by anti- clockwise, stretcher will be down

Bright arrow symbols direct the user how to operate clearly

Central-lock castors with steering, diameter 150mm

Concave plastic chassis cover, for storage use

Rotary rollover PP protective railing is blew molding by one-shot plastic forming, 3 different positions of the side rail supply different use

PP bed board produced by blowing molding, easy for clean

Standard accessories: Mattress, I.V. pole



### SPECIFICATIONS

Model	HBD8G-MS1
Back-rest	0 – 89°
Height range (bed board to ground)	791 - 509 mm
Standard accessory	Mattress, IV pole
Bed Board	1970x685 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	1970x685 mm

Net Weight	67 kg
Dimension	1960Lx680W mm

## HBD8H-CS CONNECTING TRANSFER STRETCHER

It is a manual Connecting Transfer Stretcher, can be used to defend cross-infection  
 Cart A inside the operation room, Cart B outside, adjusting them to the same height,  
 then push Cart B to Cart A, two carts will be connected, following the orbit, bedboard  
 C can be pushed from one cart to another, and VICE VERSA  
 Bright arrow symbols direct the users how to operate clearly  
 Standard accessories: Mattress, I.V. pole



### SPECIFICATIONS

Model	HBD8H-CS
Back-rest	0 – 85°
Height range (bed board to ground)	660 - 910 mm
Standard accessory	Mattress, IV pole
Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	3870x840 mm
Net Weight	139 kg

## HBD8I-HS HYDRAULIC TRANSFER STRETCHER

It is a type of hydraulic transfer stretcher  
 Five functions available by pressing the foot-pedals  
 Bright arrow symbols direct the users how to operate clearly  
 Central-lock castors, diameter 150mm. Matched with direction castor under the  
 chassis, which help users easily control the stretcher direction  
 Standard accessories: Mattress, I.V.Pole



### SPECIFICATIONS

Model	HBD8I-HS
Back-rest	0 – 85°
Height range (bed board to ground)	813 – 590 mm
Standard accessory	Mattress, IV pole

Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Trendelenburg	0 – 12°
Reverse trendelenburg	0 – 12°
Overall dimension	1906x743 mm
Net Weight	82 kg

## HBD8J-HS1 HYDRAULIC TRANSFER STRETCHER

It is a type of hydraulic transfer stretcher, 5 functions available by hydraulic and gas spring and crank system together

At the head-end, there is a pair of push- handles, can be folded up and down

Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily, Bed board for X-Ray Examination



### SPECIFICATIONS

Model	HBD8J-HS1
Back-rest	0 – 85°
Height range (bed board to ground)	860 – 630 mm
Knee Gatch	0°- 50°
Trendelenburg	0 – 13°
Reverse trendelenburg	0 – 13°
Standard accessory	IV pole
Bed Board	1932x630 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	2050x745 mm
Net Weight	105 kg
Height adjustment	835-600 mm

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)