

WHEELCHAIR



WHEELCHAIR

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of controller and other functions reflect the user-friendly design, and provide great convenience for the user.

HBD91 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



SPECIFICATIONS

Model	HBD91-MW1	HBD91-MW2	HBD91-MW3	HBD91-MW4
Wheelchair length	970 mm	1020 mm	950 mm	855 mm
Opened wheelchair width	660 mm	630 mm	540 mm	535 mm
Wheelchair Height	920 mm	890 mm	870 mm	855 mm
Closed wheelchair width	350 mm	345 mm	280 mm	250 mm
Seat Width	430 mm		405 mm	450 mm
Front wheel	6" inch			
Rear wheel	24" inch		16" inch	14" inch
Seat Height	470 mm	460 mm	445 mm	480 mm
Seat depth	400 mm		420 mm	430 mm
Backrest Height	450 mm		380 mm	350 mm
Armrest Height	210 mm		220 mm	190 mm
Distance between armrest	425 mm	445 mm	435 mm	440 mm
Height between footrests and ground	120 mm	80 mm	95 mm	60 mm
Net Weight	12.94 Kg	13.28 Kg	10.24 Kg	13.75 Kg

Model	HBD91-MW5	HBD91-MW6	HBD91-MW7	HBD91-MW8
Wheelchair length	1250 mm	1055 mm	830 mm	1180 mm
Opened wheelchair width	690 mm	650 mm	580 mm	660 mm
Wheelchair Height	1210 mm	920 mm	780 mm	1320 mm
Closed wheelchair width	365 mm	350 mm		360 mm
Seat Width	450 mm	455 mm	360 mm	460 mm
Front wheel	8" inch		4" inch	8" inch
Rear wheel	24" inch			

Seat Height	520 mm	510 mm	470 mm	520 mm
Seat depth	455 mm	450 mm	380 mm	460 mm
Backrest Height	450 mm	390 mm	330 mm	490 mm
Armrest Height	240 mm	200 mm	215 mm	230 mm
Distance between armrest	450 mm		415 mm	450 mm
Height between footrests and ground	100 mm			115 mm
Net Weight	25.4 Kg	17.4 Kg	13 Kg	23.35 Kg



HBD91-MW1



HBD91-MW2



HBD91-MW3



HBD91-MW4



HBD91-MW5



HBD91-MW6



HBD91-MW7



HBD91-MW8

HBD92 ELECTRIC WHEELCHAIR

It is a type of steady, simple and Practical Electric Wheelchair



SPECIFICATIONS

Model	HBD92-EW0	HBD92-EW1	HBD92-EW2
Wheelchair length	1055 mm	1050 mm	-
Opened wheelchair width	670 mm	685 mm	-
Wheelchair Height	880 mm	940 mm	-
Closed wheelchair width	500 mm	420 mm	-
Seat Width	445 mm		460 mm

Front wheel	8" inch		-
Rear wheel	22" inch	12" inch	-
Seat Height	480 mm	495 mm	460 mm
Seat depth	430 mm	420 mm	-
Backrest Height	400 mm	430 mm	-
Armrest Height	210 mm	205 mm	-
Distance between armrest	440 mm		-
Height between footrests and ground	100 mm		-
Km of Endurance	> 20 Km		
Gradeability (°)	> 6	-	-
Speed per Hour	> 6 Km/h		
Net Weight (Without battery)	29.1 Kg	34 Kg	-
Bearing Load	-	100 Kg	-
Gradeability	-	> 6°	-
Total Length	-	-	1000 mm
Total Width	-	-	620 mm
Total Height	-	-	890 mm
Seat Depth	-	-	400 mm
Maximum Load	-	-	160 Kg
Maximum Speed	-	-	7.2 Km/h
Mileage	-	-	17 Km
Damping System	-	-	2
Motor Power	-	-	24/250x2 V/W
Charger	-	-	24 (V2/A)
Battery	-	-	24/10 (V/Ah)
Charging Time	-	-	3-5 (h)
Front Wheel	-	-	8" inch
Rear Wheel	-	-	12" inch
Wheel Base	-	-	650 mm
Maximum Gradeability	-	-	12 °
Net Weight	-	-	25 Kg



HBD92-EW0



HBD92-EW1



HBD92-EW2

HBD93-123 ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

Thickened aluminium alloy body

strong bearing capacity

180° opening and closing, easy to realize the transfer of various scenes

Electric motor, one-button control, easy lifting operation, convenient and quick displacement

Soft and comfortable armrest, comfortable to use

360-degree universal wheel, with locking mechanism, shifting is more secure

PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract



SPECIFICATIONS

Model	HBD93-123
Material	Iron, Plastic
Maximum loading weight	150 kg
Lifting Range	25 cm, from 40 cm to 65 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 55 cm
It suites for bed	Height of bed from 11 cm to 63 cm
Waterproof level	IP44
Dimensions	123x72.5x54.5 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, rechargeable

HBD93-131 ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

Thickened aluminium alloy body, strong bearing capacity

180° opening and closing, easy to realize the transfer of various scenes

Electric motor, one-button control, easy lifting operation, convenient and quick displacement

Soft and comfortable armrest, comfortable to use. 360-degree universal wheel, with locking mechanism, shifting is more secure

PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract



SPECIFICATIONS

Model	HBD93-131
Material	Iron, Plastic
Maximum loading weight	150 kg

Lifting Range	33 cm, from 40 cm to 73 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 55 cm
It suites for bed	Height of bed from 11 cm to 72 cm
Waterproof level	IP44
Dimensions	131x72.5x54.5 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, rechargeable

HBD93-123A ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

This is a kind of patient lift
it is electric lift
height of sling can be adjustable



SPECIFICATIONS

Model	HBD93-123A
Material	Iron, Plastic, Nylon
Maximum loading weight	150 kg
Lifting Range	25 cm, from 40 cm to 65 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 55 cm
It suites for bed	Height of bed from 11 cm to 63 cm
Waterproof level	IP44
Dimensions	123x72.5x54.5 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, rechargeable

HBD93 ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

Thickened aluminium alloy body, strong bearing capacity

180° opening and closing, easy to realize the transfer of various scenes

Electric motor, one-button button control, easy lifting operation, convenient and quick displacement

Soft and comfortable armrest, comfortable to use

360-degree universal wheel, with locking mechanism, shifting is more secure

PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract



SPECIFICATIONS

Model	HBD93-198	HBD93-58
Material	Iron, Plastic, Nylon	Steel ,PU, Silicone cushion
Maximum loading weight	150 kg	
Lifting Range	Seat: 33 cm, from 40 cm to 73 cm. Upper lifting: 29 cm	Seat: 30 cm, from 40 cm to 70 cm.
Wheel	Two front wheels are with brake	
Door width, chair can pass it	At least 55 cm	At least 58 cm
It suites for bed	Height of bed from 11 cm to 72 cm	Height of bed from 9 cm to 70 cm
Waterproof level	IP44	Yes, IP44, top grade waterproof
Dimensions	198x72.5x54.5 cm	65x57x73 cm
Rated power	96 W	
Voltage	DC 24 V	
Battery	Two Batteries	-
With upper lifting vest	Yes	-
With individual soft cushion	Yes	-
Power supply	Battery, rechargeable	Battery, there is a button on battery to power it on /off
Control	-	Remote control
Handles	-	All foldable



HBD93-198



HBD93-58

HBD93-55 ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

Thickened aluminum alloy body, strong bearing capacity

180° opening and closing, easy to realize the transfer of various scenes

Electric motor, one-button control, easy lifting operation, convenient and quick displacement

Soft and comfortable armrest, comfortable to use

360-degree universal wheel, with locking mechanism, shifting is more secure

PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract



SPECIFICATIONS

Model	HBD93-55
Material	Steel ,PU
Maximum loading weight	150 kg
Lifting Range	Seat:30 cm, from 40 cm to70 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 58 cm
It suites for bed	Height of bed from 9 cm to 70 cm
Waterproof level	Yes,IP44, top grade waterproof
Control	Remote control
Handles	All foldable
Dimensions	55x65x73 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, there is a button on battery to power it on /off

HBD94 ELECTRIC TOILET LIFT

Adjustable Lifting:Adjustable lifting angle for seat height comfort at different heights
 Strong Thrust:Double push strength up to 200kg, satisfying 150kg weight and sitting freely

Safe Voltage :DC 24V safe voltage, ultra-low energy consumption

Silent Operation: The whole process of silent operation ensures the comfort of sitting

Self-protection System: When a unilateral failure occurs, the machine automatically starts the self-protection function and stops running

One-button Reset:After use, long press the switch key to reset with one button, while locking other buttons to ensure safe use



SPECIFICATIONS

Model	HBD94-150	HBD94-150A	HBD94-150B	HBD94-150C
Material	Iron, Plastic			
Seat height	From 45 to 75 cm			
Weight Capacity	150 kg			
Noise	Almost silent electric motor when operation			
Inclination angle	15°-16°			
Size	57cm width, 65cm length,47cm height			
For person	Bariatric person, elderly person,patients disabled person, and pregnant women			
Rated power	96 W / 2 A			
Voltage	DC 24 V			
Power source	Electric power	Battery power,rechargeable battery		Electric power



HBD94-150



HBD94-150A



HBD94-150B



HBD94-150C

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com