

ORTHOPEDICS TRACTOR RACK



ORTHOPEDICS TRACTOR RACK

Labstac's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.

OPT14 ORTHOPEDICS TRACTOR RACK

Adopt carbon fiber material, the prespective imaging is clear, light and strong

The plug in design is convient and quick to connect with table

Equipped with a traction frame accessory car, which is convient to move and save space



SPECIFICATIONS

Model	OPT14-1420	OPT14-1700	OPT14-1550
Table Length	1420±50 mm	1700±10 mm	1550±50 mm
Table Width	765±50 mm	750±100 mm	450±50 mm
Traction device horizontal rotary		0-180°	
Traction distance		(0-200)+10%	
Traction Device height adjustment	>400 mm	-	>190 mm



OPT14-1420



OPT14-1700



OPT14-1550

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com