

PARAFFIN TRIMMER



PARAFFIN TRIMMER

Used in Medical, Research, Tissue, Biopsy, Mega and Embedding.

Also known as Wax Trimmer, Laboratory Wax Remover, Laboratory Paraffin Remover.

PTR11-230 PARAFFIN TRIMMER

Full computer automated control.


Utilizes new type of heating elements that enable fast heating, energy, efficiency and reliability.

Temperature sensors are made in USA, precise and reliable.

Automatically saves the latest operational settings.

Memory and automatic restoration function, the preset temperature can be automatically saved after device operation.

LED displays the heating status.



SPECIFICATIONS

Model	PTR11-230
Working Area	230×160 mm
Temperature Range	Continuously adjustable within 0 - 90°C
Temperature Control Precision	±1°C
Electrical Requirement	220V±10%50 Hz or 110V±10% 50/60 Hz
Power	150 W
Dimension	255W×180D×155H mm
Weight	3 kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com