

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT11 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	CBT11-4L	CBT11-12L	CBT11-22L	CBT11-30L
Capacity	4 gallon (15L)	12 gallon (45L)	22 gallon (83L)	30 gallon (114L)
Material	Double layer anti fire cold-rolled steel PP tray			
Adjustable Shelf	1 pc		3 pcs	
Door	Single / Manual			Double / Manual
External Size	430x430x560 mm	590x460x890 mm	600x460x1650 mm	1090x460x1120 mm
Package Size	550x560x725 mm	710x580x1066 mm	720x580x1815 mm	1210x680x1285 mm
Gross Weight	38 Kg	66 Kg	116 Kg	117 Kg

Model	CBT11-45L	CBT11-60L	CBT11-90L
Capacity	45 gallon (170L)	60 gallon (227L)	90 gallon (340L)
Material	Double layer anti fire cold-rolled steel PP tray		
Adjustable Shelf	2 pcs		
Door	Double / Manual		
External Size	1090x460x1650mm	860x860x1650 mm	1090x860x1650 mm
Package Size	1210x580x1815 mm	980x980x1815 mm	1210x980x1815 mm
Gross Weight	155 Kg	195 Kg	243 Kg



CBT11-4L



CBT11-12L



CBT11-22L



CBT11-30L



CBT11-45L



CBT11-60L



CBT11-90L

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com