

STRONG ACID AND ALKALI STORAGE CABINET



STRONG ACID AND ALKALI STORAGE CABINET

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT21 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable



SPECIFICATIONS

Model	CBT21-4L	CBT21-12L	CBT21-20L	CBT21-28L
Capacity	4 gallon (15L)	12 gallon (45L)	20 gallon (75L)	28 gallon (105L)
Material	pp			
Adjustable Shelf	Optional	1 pc		2 pcs
Door	Single / Manual		Double / Manual	
External Size	430x430x500 mm	590x460x890 mm	910x590x890 mm	910x600x900 mm
Package Size	515x515x685 mm	675x545x1025 mm	995x675x1035 mm	995x685x1035 mm
Gross Weight	16 Kg	30 Kg	46 Kg	48 Kg

Model	CBT21-30L	CBT21-45L	CBT21-60L	CBT21-90L
Capacity	30 gallon (114L)	45 gallon (170L)	60 gallon (227L)	90 gallon (340L)
Material	pp			
Adjustable Shelf	2 pcs			
Door	Double / Manual			
External Size	1090x460x1120mm	1090x460x1650mm	860x860x1650 mm	1090x860x1650 mm
Package Size	1175x545x1255 mm	1175x545x1785 mm	945x945x1785 mm	1173x943x1785 mm
Gross Weight	58 Kg	72 Kg	90 Kg	107 Kg



CBT21-4L



CBT21-12L



CBT21-20L



CBT21-28L



CBT21-30L



CBT21-45L



CBT21-60L



CBT21-90L

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com