

DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER



DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

Double adjustable frequencies ultrasonic cleaner uses frequency of 45 KHz /80 KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC-18 DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

Over- Voltage and over current protection

Standard basket and lid

Low level water protection

Fully microprocessor controller and programmable

Time and temperature real time display

Standard power supply AC220+ 10%, 50/60Hz: 110V+10%, 60Hz



SPECIFICATIONS

Model	ULC18-13	ULC18-13A	ULC18-27	ULC18-27A
Capacity	13 L		27 L	
Adjustable ultrasonic frequency	45 KHz /80 KHz			
Time range/ Heating temp	1-480 mins /0- 80°C			
Ultrasonic power	300 W		500 W	
Heating power	400 W		800 W	
Tank size	300x240x180 mm		500x300x180 mm	
Packing Size	840x670x730 mm		1100x690x800 mm	
Gross Weight	73 Kg		101 Kg	

Model	ULC18-27B	ULC18-27C	ULC18-27D	ULC18-27E
Capacity	27 L			
Adjustable ultrasonic frequency	45 KHz /80 KHz			
Time range/ Heating temp	1-480 mins /0- 80°C			
Ultrasonic power	600 W		700 W	
Heating power	800 W			
Tank size	500x300x180 mm			
Packing Size	1100x690x800 mm			
Gross Weight	102 Kg		104 Kg	

Model	ULC18-3	ULC18-3A	ULC18-4	ULC18-4A
Capacity	3 L		4 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	100 W			
Heating power	400 W			
Tank size	230x140x100 mm		230x140x130 mm	
Packing Size	446x250x400 mm		446x346x453 mm	
Gross Weight	7 Kg		12 Kg	

Model	ULC18-6	ULC18-6A	ULC18-8	ULC18-8A
Capacity	6 L		8 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	200 W			
Heating power	400 W			
Tank size	300x150x150 mm		300x150x180 mm	
Packing Size	445x250x430 mm		500x350x480 mm	
Gross Weight	9.5 Kg		15.5 Kg	

Model	ULC18-10	ULC18-10A	ULC18-13B	ULC18-13C
Capacity	10 L		13 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	300 W			
Heating power	400 W			
Tank size	300x240x150 mm		300x240x180 mm	
Packing Size	446x346x483 mm		650x540x630 mm	
Gross Weight	12.5 Kg		37 Kg	



ULC18-13



ULC18-13A



ULC18-27



ULC18-27A



ULC18-27B



ULC18-27C



ULC18-27D



ULC18-27E



ULC18-3



ULC18-3A



ULC18-4



ULC18-4A



ULC18-6



ULC18-6A



ULC18-8



ULC18-8A



ULC18-10



ULC18-10A



ULC18-13B



ULC18-13C

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com