

# ULTRASONIC CLEANER



# ULTRASONIC CLEANER

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC-19 ULTRASONIC CLEANER

Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion

LED display



### SPECIFICATIONS

| Model                    | ULC19-65                 | ULC19-10       | ULC19-15       | ULC19-225      |
|--------------------------|--------------------------|----------------|----------------|----------------|
| Capacity                 | 6.5 L                    | 10 L           | 15 L           | 22.5 L         |
| Ultrasonic frequency     | Single frequency, 40 KHz |                |                |                |
| Time range/ Heating temp | 1-99 min/RT-80°C         |                |                |                |
| Ultrasonic power         | 180 W                    | 240 W          | 360 W          | 480 W          |
| Heating power            | 500 W                    |                |                |                |
| Tank size                | 300x150x150 mm           | 300x240x150 mm | 330x300x150 mm | 500x300x150 mm |
| External Size            | 330x180x310 mm           | 330x270x310 mm | 360x330x310 mm | 530x330x310 mm |
| Packing Size             | 430x250x350 mm           | 430x340x350 mm | 460x400x400 mm | 640x400x400 mm |
| Gross weight             | 6 Kg                     | 8 Kg           | 10 Kg          | 13 Kg          |
| Power Supply             | AC 220V, 50Hz            |                |                |                |

| Model                    | ULC19-30                 | ULC19-65A        | ULC19-10A      | ULC19-15A      |
|--------------------------|--------------------------|------------------|----------------|----------------|
| Capacity                 | 30 L                     | 6.5 L            | 10 L           | 15 L           |
| Ultrasonic frequency     | Single frequency, 40 KHz |                  |                |                |
| Time range/ Heating temp | 1-99 min/RT-80°C         | 1-60 min/RT-80°C |                |                |
| Ultrasonic power         | 600 W                    | 180 W            | 240 W          | 360 W          |
| Heating power            | 500 W                    |                  |                |                |
| Tank size                | 500x300x200 mm           | 300x150x150 mm   | 300x240x150 mm | 330x300x150 mm |
| External Size            | 530x330x350 mm           | 330x180x310 mm   | 330x270x310 mm | 360x330x310 mm |
| Packing Size             | 640x400x400 mm           | 430x250x350 mm   | 430x340x350 mm | 460x400x400 mm |
| Gross weight             | 15 Kg                    | 6 Kg             | 8 Kg           | 10 Kg          |
| Power Supply             | AC 220V, 50Hz            |                  |                |                |

| Model                    | ULC19-225A               | ULC19-30A      | ULC19-65B                  | ULC19-10B      |
|--------------------------|--------------------------|----------------|----------------------------|----------------|
| Capacity                 | 22.5 L                   | 30 L           | 6.5 L                      | 10 L           |
| Ultrasonic frequency     | Single frequency, 40 KHz |                | Double frequency 28/40 KHz |                |
| Time range/ Heating temp | 1-60 min/RT-80°C         |                | 1-99 min/RT-80°C           |                |
| Ultrasonic power         | 480 W                    | 600 W          | 180 W                      | 240 W          |
| Heating power            | 500 W                    |                |                            |                |
| Tank size                | 500x300x150 mm           | 300x150x150 mm |                            | 300x240x150 mm |
| External Size            | 530x330x310 mm           | 530x330x350 mm | 330x180x310 mm             | 330x270x310 mm |
| Packing Size             | 640x400x400 mm           |                | -                          | -              |
| Gross weight             | 13 Kg                    | 15 Kg          | -                          | -              |
| Power Supply             | AC 220V, 50Hz            |                | -                          | -              |

| Model                    | ULC19-15B                  | ULC19-225B     | ULC19-30B      | ULC19-65C                |
|--------------------------|----------------------------|----------------|----------------|--------------------------|
| Capacity                 | 15 L                       | 22.5 L         | 30 L           | 6.5 L                    |
| Ultrasonic frequency     | Double frequency 28/40 KHz |                |                | Single frequency, 40 KHz |
| Time range/ Heating temp | 1-99 min/RT-80°C           |                |                |                          |
| Ultrasonic power         | 360 W                      | 480 W          | 600 W          | 180 W                    |
| Heating power            | 500 W                      |                |                |                          |
| Tank size                | 330x300x150 mm             | 500x300x150 mm | 500x300x200 mm | 300x150x150 mm           |
| External Size            | 360x330x310 mm             | 530x330x310 mm | 530x330x350 mm | 300x180x310 mm           |
| Packing Size             | -                          | -              | -              | -                        |
| Gross weight             | -                          | -              | -              | -                        |
| Power Supply             | -                          | -              | -              | -                        |

| Model                    | ULC19-10C                | ULC19-15C      | ULC19-225C     | ULC19-30C      |
|--------------------------|--------------------------|----------------|----------------|----------------|
| Capacity                 | 10 L                     | 15 L           | 22.5 L         | 30 L           |
| Ultrasonic frequency     | Single frequency, 40 KHz |                |                |                |
| Time range/ Heating temp | 1-99 min/RT-80°C         |                |                |                |
| Ultrasonic power         | 240 W                    | 360 W          | 480 W          | 600 W          |
| Heating power            | 500 W                    |                |                |                |
| Tank size                | 300x240x150 mm           | 330x300x150 mm | 500x300x150 mm | 500x300x200 mm |
| External Size            | 300x270x310 mm           | 360x330x310 mm | 530x330x310 mm | 530x330x350 mm |
| Packing Size             | -                        | -              | -              | -              |
| Gross weight             | -                        | -              | -              | -              |
| Power Supply             | -                        | -              | -              | -              |



ULC19-65



ULC19-10



ULC19-15



ULC19-225



ULC19-30



ULC19-65A



ULC19-10A



ULC19-15A



ULC19-225A



ULC19-30A



ULC19-65B




ULC19-10B



ULC19-15B



ULC19-225B



ULC19-30B



ULC19-65C



ULC19-10C



ULC19-15C



ULC19-225C



ULC19-30C

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)